

自動血圧計試験結果に関する調査集計表

株式会社エー・アンド・デイ(A&D)

(医用)

製造会社名	製品名	型番	家庭用/医療用	試験結果に関する情報(学術論文等)	備考
株式会社 エー・ アンド・ デイ (A&D)	上腕型自動血圧計	TM-2430	医用血圧計	<p>Authors: Palatini P, Frigo G, Bertolo O, Roman E, Da Cortà R, Winnicki M. Article: Validation of the A&D TM-2430 device for ambulatory blood pressure monitoring and evaluation of performance according to subjects' characteristics. Published on: Blood Pressure Monitoring, 1998, 3: 255-260 Protocol: BHS Result: A/A</p> <p>Authors: Yip GW, So HK, Li AM, Tomlinson B, Wong SN, Sung RY. Article: Validation of A&D TM-2430 upper-arm blood pressure monitor for ambulatory blood pressure monitoring in children and adolescents, according to the British Hypertension Society protocol. Published on: Blood Pressure Monitoring, 2012, 17:76-9 Protocol: BHS Result: A/B</p>	(日本発売版型式: TM-2433<現行>)
株式会社 エー・ アンド・ デイ (A&D)	上腕型全自動血圧計	TM-2655P	医用血圧計	<p>Authors: Zhanna D. Kobalava, Yulia V. Kotovskaya, Lala A. Babaeva and Valentin S. Moiseev. Article: Validation of TM-2655 oscillometric device for blood pressure measurement. Published on: Blood Pressure Monitoring, 2006, 11:87-90 Protocol: BHS Result: A/A</p>	日本販売終了
株式会社 エー・ アンド・ デイ (A&D)	上腕型手動式デジタル血圧計 <スワンハート>	UM-101	医用血圧計	<p>Authors: George S. Stergiou, Periklis P. Gionas, Charilaos P. Gkinos and Dimitris G. Tzamouranis Article: Validation of the A&D UM-101 professional hybrid device for office blood pressure measurement according to the European Hypertension Society International Protocol. Published on: Blood Pressure Monitoring, 2008, 13:37-42 Protocol: ESH Result: Pass</p> <p>Authors: Pruijm MT, Wuerzner G, Glatz N, Alwan H, Ponte B, Ackermann D, Burnier M, Bochud M. Article: A new technique for simultaneous validation of two manual nonmercury auscultatory sphygmomanometers (A&D UM-101 and Accoson Greenlight 300) based on the International protocol. Published on: Blood Pressure Monitoring, 2010, 15:322-5 Protocol: ESH Result: Pass</p>	日本販売終了
株式会社 エー・ アンド・ デイ (A&D)	上腕型全自動血圧計	TM-2656	医用血圧計	<p>Authors: Wei-Fang Zeng, Ming Liu, Yuan-Yuan Kang, Yan Li and Ji-Guang Wang Article: Validation of the fully automated A&D TM-2656 blood pressure monitor according to the British Hypertension Society Protocol Published on: Blood Pressure Monitoring, 2013, 18:223-226 Protocol: BHS Result: A/A</p>	(日本発売版型式: TM-2656VPW<現行>)