

| 製造会社名 | 製品名 | 型番 | 家庭用/医療用 | 試験結果に関する情報(学術論文等) | 備考 |
|------------------------------------|-----------------------|--|---------|--|--|
| 株式会社 エー・ アンド・ デイ (A&D) | 上腕型血圧計 | UA-767 | 家庭用血圧計 | Authors: Rogoza AN, Pavlova TS, Sergeeva MV. Article: Validation of A&D UA-767 device for the self-measurement of blood pressure. Published on: Blood Pressure Monitoring, 2000, 5:227-231 Protocol: BHS Result: A/A | 日本販売終了 |
| 株式会社 エー・ アンド・ デイ (A&D) | 上腕型血圧計 | UA-767 Plus | 家庭用血圧計 | Authors: Kobalava ZD, Kotovskaia LV, Rusakova OS, Babaeva LA. Article: Validation of UA-767 Plus device for self-measurement of blood pressure. Published on: Clinical Pharmacology and Therapy, 2003, 12 (2) Protocol: BHS Result: A/A | (日本発売版型式: UA-767PC<現行>) |
| 株式会社 エー・ アンド・ デイ (A&D) | 上腕型血圧計 | UA-774 (UA-767Plus) | 家庭用血圧計 | Authors: Verdecchia P, Angeli F, Poeta F, Reboldi GP, Borgioni C, Pittavini L, Porcellat C. Article: Validation of the A&D UA-774(UA-767Plus) device for self-measurement of blood pressure. Published on: Blood Pressure Monitoring, 2004, 9:225-229 Protocol: BHS Result: A/A | (日本発売版型式: UA-774J<現行>) |
| 株式会社 エー・ アンド・ デイ (A&D) | 上腕型血圧計 | UA-787 | 家庭用血圧計 | Authors: Longo D, Toffanin G, Garbelotto R, Zaetta V, Businaro L, Palatini P. Article: Performance of the UA-787 oscillometric blood pressure monitor according to the European Society of Hypertension protocol. Published on: Blood Pressure Monitoring, 2003, 8: 91-95 Protocol: ESH Result: Pass | |
| 株式会社 エー・ アンド・ デイ (A&D) | 上腕型血圧計 | UA-85x | 家庭用血圧計 | Authors: Bonso E, Ragazzo F, Palatini P. Article: Validation of A&D UA-85x(UA-851, 852, 853, 854, 855) device for blood pressure measurement. Published on: Blood Pressure Monitoring, 2008, 13:339-341 Protocol: ESH Result: Pass | (日本発売版型式: UA-851<現行>) |
| 株式会社 エー・ アンド・ デイ (A&D) | 上腕型 手動式デジタル血圧 計 | UA-704 | 家庭用血圧計 | Authors: Zhanna D. Kobalava, Yulia V. Kotovskaya and Eugeny Rodionov. Article: Validation of semi-automatic device UA-704 for self-measurement of blood pressure. Published on: Blood Pressure Monitoring, 2005, 10:223-225 Protocol: BHS Result: A/A | |
| 株式会社 エー・ アンド・ デイ (A&D) | 上腕型血圧計 | UA-1020 (UA-1010, UA-1020, UA-1030T) | 家庭用血圧計 | Authors: Wei-Fang Zeng, Yuan-Yuan Kang, Ming Liu, Yan Li and Ji-Guang Wang Article: Validation of the A&D UA-1020 upper-arm blood pressure monitor for home blood pressure monitoring according to the British Hypertension Society Protocol Published on: Blood Pressure Monitoring, 2013, 18:177-181 Protocol: BHS Result: A/A | (日本発売版型式: UA-1020P, UA-1020G<以下、現行>, UA-1020T<販売終了>) |
| 株式会社 エー・ アンド・ デイ (A&D) | 上腕型血圧計 | UA-651 | 家庭用血圧計 | Authors: Elisabetta Benetti, Claudio Fania and Paolo Palatini Article: Validation of the A&D BP UA-651 device for home blood pressure measurement according to the European Society of Hypertension International Protocol revision 2010 Published on: Blood Pressure Monitoring, 2014, 19:50-53 Protocol: ESH Result: Pass | (日本発売版型式: UA-651BLE<現行>) |
| 株式会社 エー・ アンド・ デイ (A&D) | 上腕型血圧計 | UA-651 (with a wide-range cuff) | 家庭用血圧計 | Authors: Elisabetta Benetti, Claudio Fania and Paolo Palatini Article: Validation of the A&D BP UA-651 device with a wide-range cuff for home blood pressure measurement according to the European Society of Hypertension International Protocol revision 2010 Published on: Blood Pressure Monitoring, 2015 20:164-167 Protocol: ESH Result: Pass | (日本発売版型式: UA-651BLE<現行>) |